

NIKON FRANCE S.A.

191, Rue du Marché Rollay - 94504 Champigny sur Marne Cedex - France
 Tél. : +33-(1)-45-16-45-16 - Fax : +33-(1)-45-16-00-33 - http://www.nikon.fr


S.A. au capital de 3 820 000 Euros - RCS Créteil B 337554968 - APE 514 S

PROCEDURE D'ETALONNAGE


ISO 9001 Certified




Certificate No. 5257
 NIKON CORPORATION
 Instruments Division
 Yokohama Plant




② Eccentricity of screen center




① Magnification accuracy of projection lens




③ Parallel positioning accuracy of stage




④ Squareness of X axis and Y axis




⑦ Traveling accuracy in X axis direction



⑤ Traveling accuracy in Y axis direction



⑥ Image shift at the focus positions between while rotating Z axis handle in arrow R and L direction



⑧ Z axis accuracy

"d" consists of the errors caused by stage glass inclination and Z axis inclination and up-down guide slack.

